
Feelvita Acticur

EMS-Durchblutungsstimulator / EMS blood circulation stimulator

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GB Instruction Manual
from page 24

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+49 (0) 6431 9 21 96 00

Content

Product components	24
General information	25
Explanation of symbols	25
Safety.....	25
Proper use	25
Safety instructions	25
Precautions during use	26
Before first use	28
Using the battery	29
Remote control battery	29
Machine battery.....	29
Application.....	29
EMS mode	30
TENS mode	32
Cleaning and care.....	33
Troubleshooting	34
EMC Statements	36
Symbols	36
Executive Standards	38
Appendix	38
Disposal	41
Technical data	42

Product components

- 1 Main unit
- 2 Foot plates
- 3 On/off button
- 4 Charging socket
- 5 Connection for electrode cables
- 6 FlexShifter
- 7 Carrying handle
- 8 Receiver for remote control signal
- 9 Remote control
- 10 AAA battery (2x)
- 11 Electrode pad, small (2x)
- 12 Electrode pad, large (1x)
- 13 Electrode cable
- 14 Moisturising cream
- 15 Power supply
- 16 Charging cable

PLEASE TAKE THE TIME TO READ THESE INSTRUCTIONS FOR USE CAREFULLY.

FOR YOUR OWN PROTECTION AND THE PROTECTION OF THIRD PARTIES, PAY PARTICULAR ATTENTION TO THE CHAPTER ON SAFETY INSTRUCTIONS AND THE WARNINGS CONTAINED IN THE DOCUMENT.

KEEP THESE INSTRUCTIONS FOR USE FOR FUTURE REFERENCE. ONLY PASS THE DEVICE ON TO THIRD PARTIES TOGETHER WITH THE OPERATING INSTRUCTIONS.

THE ILLUSTRATIONS IN THESE INSTRUCTIONS FOR USE ARE INTENDED AS A SCHEMATIC REPRESENTATION OF THE CORRECT APPLICATION. THEY MAY DIFFER IN DETAIL FROM THE ORIGINAL PRODUCT.

All available accessories for the appliance are included in these operating instructions. The scope of delivery of the set you have purchased can be found on the packaging or on the separate insert. Parts not included in your set can of course be purchased at a later date.

General information

Explanation of symbols

The following symbols and signal words are used in this manual, on the device or on the packaging.

WARNING!

Failure to observe this warning could lead to a serious or potentially lethal accident.

CAUTION!

Failure to observe this caution could lead to injury.

NOTE!

Information about potential material damage.

Safety

Proper use

- The EMS blood circulation stimulator is designed to relax and regenerate muscles using electrical impulses and to promote blood circulation. It is primarily intended for use on the feet and legs. Other areas of application are possible when using original accessories. It is essential to observe the restrictions on use in the „Safety instructions“ section!
- The device is intended for domestic use only. It is not suitable for commercial or industrial use.
- Only use the device as described in the instructions. Any other use is considered improper.

Safety instructions

- This device must not be used by children. Children must not play with the device.
- The device is not suitable for use by persons with reduced physical, sensory or mental capabilities, unless they are supervised by a person responsible for their safety.

WARNING!

RESTRICTIONS ON USE – to avoid health risks, the following persons must NOT use the device!

- Patients with local pacemakers, metal implants or an automatic implantable cardioverter defibrillator (ALCD)
- Women during pregnancy
- Patients with local acute purulent inflammation
- Patients with local malignant tumours

- Patients with skin sensitivity disorders
- Patients with suspected or diagnosed heart disease
- Patients who have symptoms of existing deep vein thrombosis (DVT) or are being treated for this form of thrombosis

Thrombosis can also develop without any symptoms. However, it is important to recognise the typical signs:

- Swelling, pain and tenderness in one of the legs (calf)
- Severe pain in the affected area
- Redness of the skin on the back of the leg below the knee
- Warm skin at the site of the blood clot

If these symptoms occur, seek medical attention immediately.

Consult your doctor before use if:

- You are currently under medical treatment
- You currently have or have had heart problems
- You are currently undergoing pain treatment
- you suffer from a tendency to bleed
- You suffer from suspected or diagnosed epilepsy
- You are unsure of the cause of your symptoms or whether the EMS blood circulation stimulator is suitable for you

Do not use the device:

- In the presence of electronic monitoring devices
- In conjunction with a life-supporting electronic device
- In the bathroom or in connection with water
- While you are sleeping
- In connection with a splint or plaster cast without first consulting a doctor

Precautions during use

Foot surfaces and electrode pads

- DO NOT use the electrode pads on the head, face, neck or chest area
- Placing the electrode pads near the chest may increase the risk of cardiac arrhythmia
- Only apply to areas with intact/undamaged skin
- DO NOT apply electrode pads on or near the following areas:

- Cuts and injuries
 - Open wounds or areas with skin rash
 - Swollen, reddened, infected or inflamed areas
 - Skin eruptions, e.g. phlebitis, thrombophlebitis, varicose veins or cellulitis
 - On the genitals
 - On areas that have been treated with radiotherapy (within the last 6 months)
 - On opposite thighs (see also instructions for use)
 - Over the carotid sinus nerve (front of the neck)
- Do not use the electrode pads on areas of skin that are insensitive. Skin irritation may occur if the intensity is too high.
 - The electrode pads are for use by one person only. Do not share with other people.
 - Skin irritation may occur under the electrode pads.
 - During the initial phase, symptoms may worsen at first, as the body needs to get used to the increased muscle stimulation and blood flow. In this case, reduce the intensity and duration of the treatment.
 - If you experience any unwanted side effects or are unsure about the application or symptoms, please consult your doctor.

**WARNING!****RISK OF ELECTRIC SHOCK OR FIRE**

- Check the power cord for damage before use. Do not use damaged power adapters or power cords. Only connect the power adapter to sockets that carry the electrical voltage specified on the power adapter.
- If the power cord of the device is damaged, it must be replaced with a special connection cable available from the manufacturer or its customer service department.
- Do not lay the power cord over sharp metal edges, do not place heavy objects on it and protect it from heat. Do not kink, twist or knot the power cord.
- Never immerse the device, the mains plug or the mains cable in water or other liquids. Do not touch the device, the mains cable or the mains plug with wet hands.
- If the device falls into water, unplug it immediately. Never attempt to pull the device out of the water while it is connected to the mains.
- Only charge the battery with the supplied power adapter. Do not charge any other devices with the supplied power adapter.
- The device cannot be used while it is connected to the power supply unit or the mains.

- Keep the device away from heat sources (e.g. hotplates or ovens) and open flames to prevent damage.
- Only use and store the device in dry indoor areas.
- In case of fire, do not extinguish with water! Smother flames with a fire blanket or a suitable fire extinguisher.

Danger from rechargeable batteries/rechargeable accumulators:

- Improper handling of batteries and rechargeable batteries may result in an explosion. Contact with leaked battery acid may cause chemical burns.
- The device has a permanently installed Li-ion battery.
- Do not heat the battery, short-circuit it or throw it into an open fire.
- Avoid contact of battery acid with skin, eyes and mucous membranes. In case of contact with the acid, rinse the affected areas immediately with plenty of clean water and seek medical attention if necessary.

NOTE

Risk of damage

Improper handling of the device may result in damage.

- Do not stand on the device. Only use it while seated.
- Do not expose the device to extreme temperatures. Protect it from direct sunlight.
- Store the device in a cool, dry place out of reach of children and animals.
- Do not use the device if the plastic parts are cracked or broken or if the device is damaged in any other way.
- Only use the device with original accessories, otherwise you will damage the device and compromise its safety.

Before first use

CAUTION!

Risk of suffocation from plastic film and bags! Keep children and animals away from the packaging material.

- **Check** the set for completeness and any transport damage. The scope of delivery of the set you have purchased can be found on the packaging or in the separate insert. If you have any complaints, please contact our customer service.
- **Remove** the packaging material and any labels on the product that are used for transport protection or advertising
- **Insert** the batteries into the remote control.

Charging the device

1. Make sure that the device is switched off.
2. Connect the charging cable to the charging socket (FlexShifter on the right-hand side) and the power supply unit.
3. Connect the power supply unit to a working power outlet.
4. During charging, the on/off button on the device flashes white.
5. When the battery is fully charged, the on/off button will light up white continuously.

Using the battery

Remote control battery

1. The remote control uses 2 AAA batteries, which could be used for about half a year (if using for 30 minutes a day).
2. When „“, is displayed on the remote control screen, it means that the battery of the remote control is low, please replace it with a new battery.
3. If the product is not used for a long time, please take out the battery, otherwise the battery may leak and damage the machine.
4. Do not use batteries other than specified.
5. Pay attention to the positive and negative polarity when installing the battery, „+“ is positive, „-“ is negative.
6. Do not use this battery when the temperature exceeds more than 45°C, otherwise it will affect the performance and life of the battery.
7. Please dispose of flat batteries according to your local environmental regulations.

Machine battery

1. When you hear a „beep“ 3 times in a row when using the machine, it means the lithium batteries are low and need to be charged using the power adapter. Use the supplied charging cable, with one end plugged into the DC 5.0V and the other end into the charging port of the machine. It takes approximately 3 hours to fully charge.
2. A fully charged battery will last for approximately 60 days of use at 30 minutes of use per day.
3. Under normal use, the service life of the rechargeable batteries is at least 500 times of charges and discharges for approximately 2-3 years (using twice per day for no more than 30 min each time). If the batteries are not used for a long time, their service life will be shortened so it's recommended to charge at least once every three months.

Application

Remote

- The remote control can be used to control the various modes/intensities.
- Point the remote control at the receiver in the recess where the remote control can be stored.
- Activate the remote control with the on/off button.
- Pressing and holding the on/off button switches the remote control and device off when a connection is established.

On/off button, also used to select the preset time

Selecting preset programmes

„Foot“ button: Increase/decrease the intensity of the foot pads

„Body“ button: Increase/decrease the intensity of the electrode pads

Timer

- When the remote control is activated, the preset programme duration can be adjusted by repeatedly pressing the on/off button (15/30/45/60 minutes). The time is shown on the display.
- The recommended treatment time is 30 minutes daily.

Modes

- By pressing the mode button , you can choose between 9 preset treatment programmes.

Intensity

- Depending on the application, the strength of the electrical pulses can be adjusted using the +/- buttons.
- „Foot“ button: Increase/decrease the intensity of the foot pads.
- Body button: Increase/decrease the intensity of the electrode pads.
- The device has 20 intensity levels:
 - 1-5 low
 - 6-10 medium
 - 11-15 high
 - 16-20 very high

EMS mode

The following should be achieved by using the foot pads:

- Reduction of pain and discomfort in the feet, legs and/or ankles caused by circulatory disorders (peripheral arterial disease)
- Improvement of blood circulation to prevent or reduce blood congestion caused by circulatory disorders (chronic venous insufficiency/varicose veins)
- Improvement of blood circulation to prevent or reduce blood congestion caused by immobility after surgery on the lower limbs or by diabetes
- Improvement of symptoms associated with chronic venous insufficiency/varicose veins
- Supports leg vein health by improving blood circulation, increasing oxygenation of the blood and reducing swelling in the feet, legs and ankles

The following can be supported by using the device if you sit a lot or are inactive for long periods of time:

- Relief of tired, heavy or painful legs, including cramps
- Support for leg vein health
- Reduction of swelling in feet and ankles
- Strengthening of the lower leg muscles
- Improvement of blood circulation

Application to the soles of the feet (EMS mode)

1. Clean your feet and the foot pads with a slightly damp towel.
2. Check the underside of your feet for injuries.
3. Sit down comfortably. Your legs should form a right angle.
4. Place the device on the floor in front of you so that you can read the logo.
5. **CAUTION** – The device is only suitable for use while seated. Never stand on the device.
6. Place both feet barefoot on the corresponding surfaces. The raised areas in the middle of the device serve as a guide for the optimal positioning of the soles of your feet.
7. To switch on the device, first press the on/off button on the main unit. You will hear a beep and the display will flash turquoise.
8. Then press the on/off button on the remote control to switch it on as well.
9. You can use the remote control to operate the device and set the various modes.
10. The device starts muscle stimulation as soon as an intensity is selected on the corresponding button (Foot).

NOTE:

Select a lower intensity during the familiarisation phase.

After a short time, you should feel the muscle contractions. Depending on the programme and intensity, this can take up to 1 minute.

If no muscle contraction occurs, increase the intensity carefully and gradually. It is normal to feel a tingling sensation or varying degrees of discomfort in your feet and lower leg muscles.

After a few days of getting used to the device, you can slowly increase the intensity. However, you should still feel a pleasant contraction in your lower leg muscles.

FlexShifter function

- The automatic muscle contractions cause your feet to want to move up and down. To allow this movement, the device has a FlexShifter function.
- This ensures that the device adapts optimally to the movement rhythm of your feet, whereby the change in the angle of the joints in the foot and ankle during the electrical stimulation cycle improves venous blood flow.
- The resulting movement of the ankles caused by the FlexShifter indicates that the intensity is sufficient for therapy. Further increases improve the therapy, but should never cause discomfort.

Skin detection function

- The device has sensors to detect the moisture content of the skin. These are located under the foot surfaces.
- The measurement begins as soon as both feet are placed on the device after switching it on and an intensity of at least level 1 has been selected.
- For dry skin:
 - The operating indicator on the device changes from turquoise to green.
 - Automatic increase in intensity

- For normal skin:
 - The display lights up turquoise
 - The device operates according to the standard programme
- After 1 minute, skin detection is switched off and the indicator lights up turquoise in all cases. If the feet are removed and repositioned, skin detection starts again from the beginning

After the preset time has elapsed, „00“ appears on the display on the remote control and a beep sounds.

To switch off, press and hold the on/off button on the remote control or press the on/off button on the main unit.

NOTE:

The main unit switches off automatically if no switch-on signal is received from the remote control. To reactivate the main unit, press the on/off button.

TENS mode

The following should be achieved by using the electrode pads:

- Temporary relief of musculoskeletal pain or neuropathic pain
- Relief of chronic pain associated with musculoskeletal pain
- Relief of arthritis-related pain
- Reduction of pain and discomfort in the feet, legs and ankles caused by diabetic peripheral neuropathy


NOTE!

The small electrode pads and the large electrode pad are for use by one person only. Do not use at the same time as another person.

WARNING!

DO NOT USE THE PADS ON THE HEAD, NECK OR CHEST AREA.

Using the electrode pads (TENS mode)

1. Clean your skin thoroughly and dry it before using the electrode pads.
2. Make sure that the device is switched off.
3. Attach the electrode pads to the electrode cable and connect the other end of the electrode cable to the designated opening on the main unit (FlexShifter on the right side).
4. Carefully peel the electrode pads off the plastic film and stick the pads to the area of the body to be treated.

Neck

Shoulder

Back

Calf

Hand

Hip

Joint

Position the pads so that the area of greatest pain is between the electrode pads.
Never place the electrode pads directly on a joint.

1. To switch on the device, first press the on/off button on the main unit. You will hear a beep and the display will flash turquoise.
2. Then press the on/off button on the remote control to switch it on as well.
3. You can use the remote control to operate the device and set the different modes.
4. The device starts muscle stimulation as soon as an intensity is selected on the corresponding button (Body).
5. After the treatment is complete, remove the electrode pads from your skin. Do not pull on the electrode cable, as this could damage it.
6. After each use, stick the electrode pads back onto the plastic film.
7. Disconnect the electrode cable from the side socket of the device and store it together with the other accessories.
8. The electrode pads should always be kept clean and free of dirt, oil, sticky substances, etc., as otherwise their adhesive properties will be impaired.

NOTE:

The electrode pads can be used at the same time as the foot pads. The intensity can be controlled separately using the corresponding button on the remote control.

Select a lower intensity during the initial phase.

After a short time, you should feel muscle contractions. Depending on the programme and intensity, this can take up to 1 minute.

If no muscle contractions occur, increase the intensity carefully and gradually. It is normal to feel a tingling sensation or varying degrees of discomfort.

After a few days of getting used to the device, you can slowly increase the intensity. However, you should still feel a pleasant contraction in the respective body part.

Cleaning and care

WARNING!

Risk of electric shock

Never hold the device under running water and do not immerse it in water or other liquids.

Device / foot pads:

- Clean the device / foot surfaces by wiping them with a damp cloth.
- Do not use sharp or abrasive cleaning agents or cleaning pads.
- Store the device in a dry, cool place away from direct sunlight and out of the reach of children and animals.

Electrode pads:

- After use, stick them back onto the plastic film provided.
- If necessary, clean under running water after use and allow to dry before storing/reusing.
- Last for up to 100 applications; replacement pads can be purchased at www.genius.tv.
- Do not disconnect from the electrode cable for storage.

Ambient conditions

(Transportation and storage environment between secondary uses)

	Temperature	Humidity	Barometric pressure
Normal operation	+5°C - +40°C	15%- 93% RH	700hPa - 1060hPa
Storage	-20°C - +55°C	00%- 93% RH	700hPa - 1060hPa
Transport	-20°C - +55°C	15%- 93% RH	700hPa - 1060hPa

When the ambient temperature is 20 °C, the time required for the product to be ready for its intended use from the minimum storage temperature after use is about 2 hours; At an ambient temperature of 20 °C, the time required for the product to go from the maximum storage temperature after use to the time the product is ready for its intended use is approximately 2 hours.

Troubleshooting

Problem	Possible cause	Solution
The EMS blood circulation stimulator is switched on, but I cannot feel the electrical impulses through the foot pads.	Your bare feet are not on the foot pads at the same time.	Make sure that you are barefoot and that both feet are on the foot pads at the same time. Increase the intensity level on the remote control to the maximum level 20 so that you can feel the stimulation.
	Your soles may be too dry.	Apply a moisturising cream to the soles of your feet to improve electrical conductivity and the effect of the stimulation. Then start the application again. It may also be necessary to increase the intensity level.
	Your body may be dehydrated.	Make sure you drink enough water before and after use. The device uses the body to close the electrical circuit. Water conducts electricity very well. Dehydration increases skin resistance, which can weaken the stimulation. It is therefore essential to stay well hydrated.
	The intensity level may be too low.	The device is very safe to use. Increase the intensity gradually up to level 20 until you feel an effect. A higher level may feel unfamiliar or uncomfortable at first. The aim is not to reach level 20 immediately, but to find an intensity that is comfortable for you.

	If you do not feel any stimulation despite taking all these measures:	Test the device by placing both hands on the heel area of the foot pads at the same time. Switch on the EMS blood circulation stimulator and gradually increase the intensity to the maximum level 20 or until you feel a clear stimulation in your hands. If this is the case, the device is working correctly. If you still do not feel anything at the highest level, please contact your official distributor.
The on/off button on the EMS blood circulation stimulator remains off and the device does not work even though it is switched on.	The power supply may not be properly connected to the power outlet or the device (device does not charge).	Make sure that the power supply unit is properly plugged into the power outlet and connected to the device. If the problem persists, please contact your authorised dealer.
	The device battery is empty.	Connect the device to a power outlet to charge the battery or operate it directly from the mains.
The EMS blood circulation stimulator does not vibrate.	The device is not designed to vibrate.	
The FlexShifter function does not perform a rolling motion.	The movement is caused by muscle activity in the legs – the FlexShifter merely serves as a pivot point.	The FlexShifter only moves at higher intensities, as this activates the lower leg muscles. You may not be able to increase the intensity significantly right away – allow yourself time to get used to the sensation. It is important to select an intensity level that is comfortable for you.
	The device may be positioned too far away from you.	Make sure that your sitting position is correct: your knees should be at right angles.
The FlexShifter function makes a knocking noise on hard surfaces.	The device may not be positioned correctly or the intensity may be set too high.	Reposition the EMS blood circulation stimulator or reduce the intensity. This should make the noise quieter. Alternatively, a floor mat can be placed underneath to dampen the knocking noise.
Pain occurs in the legs after use.	The intensity level was probably too high, causing the muscles to be overworked.	Allow your muscles sufficient recovery time between applications (at least 48 hours). Start with a lower intensity and increase it slowly until the stimulation is pleasantly noticeable – without exertion. If necessary, also reduce the application time until your muscles have become accustomed to the treatment.
The stimulation with the electrode pads feels uncomfortable.	The intensity may be too high or the skin in that area may be irritated or injured.	Reduce the intensity level on the remote control using the "-" button. If your skin is irritated, wait until it has completely healed before using the device.

During use of the electrode pads, pain occurs in the thigh muscles before the session is over.	The muscles may not yet be strong enough for a full treatment.	Switch off the device early and end the session. The next time you use the device, select a lower intensity level at which you feel a pleasant, gentle stimulation. Also shorten the treatment time until your muscles get used to the stimulation.
--	--	---

CAUTION

If the problem persists, please contact our customer service. Never open the device!

EMC Statements

WARNING!

**USING ACCESSORIES OTHER THAN THOSE SPECIFIED MAY RESULT IN INCREASED EMISSIONS OR DECREASED IMMUNITY OF THE PRODUCT.
GUIDANCE AND MANUFACTURER'S DECLARATION CAN BE SEEN IN THE ATTACHMENT.**

1. The product requires special attention to EMC and requires to be installed.
Please use it based on the EMC information provided, since it may be affected by portable and mobile RF communication devices.
2. The product has been thoroughly tested and inspected to ensure its excellent performance and proper operation.

NOTE:

The product should not be used near or on any other device; if it is necessary for this, please observe if it can work normally.

Following cable types must be applied to ensure compliance with standard of radiation interference and noise immunity

Cables	Length (m)	Whether shielded
Charging cable	0.5/1	No
Electrode connecting wire	1.2	No

Symbols

Notation	Meaning
	Batch/lot number
	Serial No.
	Manufacturer

	European authorised representative
	Date of manufacture
	Warning
	Type BF applied part
	Medical device
	Warning on danger
	Category II medical device
	It means prohibition (something that must be prohibited)
	It means compulsory execution (something that must be followed)
	Indoor use only
	No sun exposure
	No rainwater
	„WEEE (Waste Electrical and Electronic Equipment)“. Wastes shall be disposed in accordance with laws.
	CE Mark and Agency Code
	Standby
	Please refer to the user manual
	Temperature limitations. It indicates the range of temperature to which medical devices can be safely exposed. The upper and lower temperature limits shall be indicated near the upper and lower horizontal lines.

	Humidity limits. It indicates the range of humidity to which medical devices can be safely exposed. The humidity limits shall be indicated near the upper and lower horizontal lines.
	Atmospheric pressure limits. It indicates the range of atmospheric pressure to which medical devices can be safely exposed. The atmospheric pressure limits shall be indicated near the upper and lower horizontal lines.
	Direct Current
	Electrodes
IP21	Dustproof and waterproof level, to avoid the falling of solid objects larger than 12.5mm and the vertical penetration of water to cause adverse effects.

Executive Standards

The product should be in accordance with the following standards and regulations:

1. IEC 60601-1:2005+A1:2012 Medical electrical equipment-Part 1: General requirements for basic safety and essential performance.
2. IEC 60601-1-11:2020 Medical Electrical Equipment - Part 1-11:
General Requirements For Basic Safety And Essential Performance Collateral Standard: Requirements For Medical Electrical Equipmentd And Medical Electrical Systems Used In The Home Healthcare Environment.
3. IEC 60601-2-10: 2012+A1:2016 Medical electrical equipment - Part 2-10: Particular requirements for the safety of nerve and muscle stimulatorss
4. IEC 60601-1-2:2014+A1:2020 Medical Electrical Equipment - Part 1-2: General Requirements For Basic Safety And Essential Performance - Collateral Standard: Electromagnetic Disturbances - Requirementsd And Tests.

Appendix

Guidance and manufacture's declaration – electromagnetic emission		
The EMS blood circulation stimulator is intended for use in the electromagnetic environment specified below. The buyer or user of the EMS blood circulation stimulator should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The EMS blood circulation stimulator use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference to nearby electronic equipment.

RF emission CISPR 11	Class B	The EMS blood circulation stimulator is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	N/A	

Guidance and manufacture's declaration – electromagnetic immunity

The EMS blood circulation stimulator is intended for use in the electromagnetic environment specified below. The buyer or user of EMS blood circulation stimulator should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+8kV contact +2kV, + 4kV. +8kV, + 15kV air	+8kV contact +2kV, + 4kV. +8kV, + 15kV air	Floors should be wood, concrete or ceramic tile. If floor is covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1kV line(s) and neutral	± 1kV line(s) and neutral	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 s	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the EMS blood circulation stimulator requires continued operation during power mains interruptions, it is recommended that the EMS blood circulation stimulator be powered from an uninterruptible power supply.

Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
--	------	------	--

NOTE:

UT is the a.c. mains voltage prior to application of the test level.

Recommended separation distances between portable and mobile RF communications equipment and the EMS blood circulation stimulator.

The EMS blood circulation stimulator is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The buyer or user of the EMS blood circulation stimulator can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the EMS blood circulation stimulator as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 KHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For rated maximum output power of transmitter that's not listed in the above chart, the recommended isolation distance d (unit is meter) can be determined by equation in the corresponding transmitter frequency column. P is the maximum output rated power of transmitter provided by its manufacturer, the unit is watt (W).

NOTE 1:

When under 80MHz and 800MHz, use equation of higher frequency.

NOTE 2:

These guidelines couldn't possibly fit all situations. Electromagnetic transmission is affected by buildings, objects, the absorption and reflection of human body.

Shenzhen Kentro Medical Electronics Co., Ltd
 2nd Floor, No. 11, Shanzhuang Road, Xikeng Village, Yuanshan Street,
 Longgang District, Shenzhen City, Guangdong Province, China
 [Tel]: +86(755)33825998, [Fax]: +86(755)33825996

Share Info GmbH
 Am Schulzentrum 12, 41564 Kaarst, Germany
 Tel: 0179 5666 508
 Email: eu-rep@share-info.com

Disposal

Packaging

The packaging material is recyclable. Dispose of the packaging in an environmentally friendly manner and take it to the recycling centre.

Electrical appliances

The crossed-out wheellie bin symbol indicates that this product must not be disposed of with normal household waste. Please dispose of electrical and electronic devices via your local take-back system. For more information on the disposal and recycling of this appliance, please contact your local waste disposal service, the shop where you purchased the product or the manufacturer.

Rechargeable batteries and batteries must be removed from the device before disposal and disposed of separately from the device. For environmental reasons, batteries and rechargeable batteries must not be disposed of with normal household waste, but must be taken to appropriate collection points. Please check your local regulations for the separate disposal of rechargeable batteries and batteries.

Remove the rechargeable battery

ATTENTION!

Do not open the battery, as there is a risk of short circuit. In addition, irritating vapours or corrosive liquids may escape.

Do not expose the battery to strong sunlight, heat or fire, as there is a risk of explosion.

Remove the battery from the device and dispose of it properly.

- To do this, remove the screw covers on the underside of the device with a pointed object.
- Loosen the screws and remove the cover.
- You can now remove the battery.
- Even discharged batteries still contain residual energy. Therefore, cover the terminals with insulating tape before disposal to prevent a short circuit, which could cause a fire or explosion.

Technical data

DEVICE	Transcutaneous Electrical Nerve Stimulator
Model	KTR-4021
Frequency range	1-100 Hz
Diameter	approx. 36.5 cm
Height	approx. 7 cm
Power	5
Battery	3.7 V; lithium battery
Battery capacity	2,200 mAh
Battery charging time	3 hours
Operating time	max. 30 hours
Remote control	2 x AAA batteries
Weight	approx. 1950 g
POWER	
Power supply	
In	100 – 240 V ~50/60 Hz
Out	5 V / 1 A
Protection	II

Pulse width	100-400 μ s
DC component	0V
Peak output voltage	5-80V
Maximum output amplitude RMS	<10V
Peak output voltage during open circuit measurement	$\leq 500V$
Output amplitude adjustment	Be able to withstand the effects of open and short circuits at output terminal, with no degrading in performance.
Adjustment of output amplitude	Continuous and uniform, keeping the minimum output not greater than 2% of the maximum output
Rated load impedance	500 Ω , with $\pm 10\%$ error
Set treatment time	15/30/45/60minutes, with an acceptable error of $\pm 10\%$

Describe for software:

Software name: KTR-4021-2081A V1.0

Version: V1.0

Customer service

If you have any questions about the device or spare parts/accessories, please contact your local retailer.

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Shenzhen Kentro Medical Electronics Co.,Ltd